



HemoWeld - BATTERY

COMPACT - LIGHT - QUICK




DELCON

Our technology. Your health.

KEY FEATURES

6



HemoWeld-Battery is the easiest way to perform sterile sealing procedures.

- 1. Portable** - supplied with a comfortable shoulder bag
- 2. Lightweight** - less than 1 kg
- 3. Lightning quick** - allows for an almost instantaneous sealing motion
- 4. Powerful** - can perform more than 500 procedures
- 5. Long-lasting** - rechargeable battery pack
- 6. Easy snap-apart separation** - a clear-cut notch line makes for an easy separation of the tube
- 7. Ergonomic handle** - with integrated splash-guard and tube retention
- 8. Anti-spark & overheat protection systems** - to avoid damage to system and tubing
- 9. Alerts notifications** - a beep and LED signals to indicate possible problems

7



9



Delcon reserves the right to modify, add, and remove the functional and visual features of the devices, at its sole discretion, at any time, and without prior notice

HemoWeld-Battery is a **compact** and **lightweight** tube sealer, easy to use and to clean, supplied together with a comfortable shoulder bag: it can be used to seal most of standard tubing available on the market. Its new built-in RF electronics **optimize the power consumption**, allowing more than 500 procedures with a fully charged battery.

Each seal by HemoWeld-Battery has a **notch** facilitating the snapping-apart of the two parts. The anti-spark and overheat **protection** system, standardly included, avoids any possible damage to the system and to the tubing. The **tube positioning system** is integrated into the **splash guard**. Any possible presence of sparks while sealing, dirty electrodes or low battery is indicated by a beep and LED signals.

HemoWeld-Battery is **certified** to work even with the needle still connected to the patient/donor.

Acoustic Alerts

Together with a dynamic LED signal, a beep indicates any problem, such as sparks during the sealing process, dirty electrodes and the low-level charge of the battery.

Ergonomic Handling

A redesigned body shape fit perfectly in your hand, making the sealing motion natural and comfortable. The new splash guard is safe for the operator and easy to remove and clean.

Long Lasting Battery

The built-in RF electronics optimize the power consumption, which means more consecutive sealings. The battery pack is rechargeable and each charge will last more than 500 procedures.

Sales Codes

HWELD-B.DB.00	Mobile tube sealer for blood bags tubing (includes RF generator, battery charger, shoulder bag and spare coaxial cable)
---------------	--

Regulation

CE marking:	Complies with Class 1 Medical Device according to directive 93/42 MDD and 2007/47 MDD
Safety Class:	1
Applied part:	CF



HEMOWELD - B

IS PART OF THE HEMOWELD LINE

A range of products specifically designed to produce sterile seals on standard PVC tubing normally used in conjunction with blood bags.

Discover all the HemoWeld Line on our website!



REV.05 - 2019

For ordering and technical information, please download the full DataSheet from our website.

DELCON srl
Via Matteucci, 25/2720862
Arcore (MB) - Italy
Phone +39 039.617670

